

DCTD Standard Operating Procedures (SOP)

Title:	γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides	Page 1 of 29
Doc. #:	SOP340543	Revision: A
		Effective Date: 03/29/2019

National Clinical Target Validation Laboratory

Applied/Developmental Research Directorate, Leidos Biomedical Research, Inc.

Frederick National Laboratory for Cancer Research

Technical Reviewer: Lindsay Dutko Date: _____

NCTVL Approval: Jiuping Ji Date: _____

IQC Approval: Katherine V. Ferry-Galow Date: _____

LHTP Approval: Ralph E. Parchment Date: _____

DCTD OD Approval: Toby Hecht Date: _____

Change History

Revision	Approval Date	Description	Originator	Approval
A	3/29/2019	Change γH2AX conjugate from FITC to biotin; change Anti-DIG-AF647 staining procedures; minor updates due to software updates.	KFG/LD	REP
--	5/25/2017	IFA staining using Bond RX Automatic Slide Staining System with γH2AX, pNBS1 and β-Catenin for tumor segmentation	DK/KFG	JJ

DCTD Standard Operating Procedures (SOP)

Title:	γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides			Page 2 of 29
Doc. #:	SOP340543	Revision:	A	Effective Date: 03/29/2019

TABLE OF CONTENTS

OVERVIEW OF IMMUNOFLUORESCENCE ASSAY FOR BIOPSIES3

1.0 PURPOSE4

2.0 SCOPE.....4

3.0 ABBREVIATIONS4

4.0 INTRODUCTION.....4

5.0 ROLES AND RESPONSIBILITIES.....5

6.0 MATERIALS AND EQUIPMENT REQUIRED6

7.0 OPERATING PROCEDURES8

8.0 OPTIONAL: SHIP TO CERTIFIED ASSAY SITE FOR ANALYSIS.....20

APPENDIX 1: BATCH RECORD21


APPENDIX 2: BOND-RX PROCESSING MODULE.....26

APPENDIX 3: STAINED SLIDES SHIPPING MANIFEST29

DCTD Standard Operating Procedures (SOP)

Title:	γ H2AX, pNBS1 IFA Staining with β -Catenin Segmentation for Tumor Biopsy Slides			Page 3 of 29
Doc. #:	SOP340543	Revision:	A	Effective Date:
				03/29/2019

OVERVIEW OF IMMUNOFLUORESCENCE ASSAY FOR BIOPSIES



DCTD Standard Operating Procedures (SOP)

Title:	γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides			Page 4 of 29
Doc. #:	SOP340543	Revision:	A	Effective Date:
				03/29/2019

1.0 PURPOSE

To standardize immunohistochemical methods for staining formalin-fixed paraffin-embedded (FFPE) tissue biopsy sections to detect and quantify histone H2AX phosphorylated at serine 139 (γH2AX) and NBS1 phosphorylated at serine 343 (pNBS1) using β-Catenin tumor segmentation for pharmacodynamic (PD) evaluation of DNA damage repair status. The goal of the SOP and associated training is to ensure consistency of biomarker measurements between operators and clinical sites.

2.0 SCOPE

This procedure applies to all personnel involved in using γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides from patients participating in clinical trials. This SOP outlines the recommended procedure for staining slides made from paraffin-embedded tumor biopsy sections using the Leica Bond-RX™ Automatic Staining System.

3.0 ABBREVIATIONS

Ab	=	Antibody
DAPI	=	4',6-Diamidino-2-Phenylindole
DCTD	=	Division of Cancer Treatment and Diagnosis
DI	=	Deionized
ER	=	Epitope Retrieval
FFPE	=	Formalin-fixed paraffin-embedded tissue
γH2AX	=	Histone H2AX Phosphorylated at Serine 139
H&E	=	Hematoxylin and Eosin
HIER	=	Heat-Induced Epitope Retrieval
ID	=	Identification/Identifier
IFA	=	Immunofluorescence Assay
LHTP	=	Laboratory of Human Toxicology & Pharmacology
NA	=	Numerical Aperture
NCTVL	=	National Clinical Target Validation Laboratory
PBS	=	Phosphate-Buffered Saline
pNBS1	=	NBS1 phosphorylated at serine 343
QC	=	Quality Control
SOP	=	Standard Operating Procedure

4.0 INTRODUCTION

The γH2AX, pNBS1 IFA with β-Catenin segmentation is an immunohistochemistry-based staining assay developed to quantify γH2AX and pNBS1 using β-Catenin staining to segment tumor from non-tumor to precisely limit quantitation of these DNA damage response markers to the tumor cells within a tumor biopsy section.

DCTD Standard Operating Procedures (SOP)

Title:	γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides			Page 5 of 29
Doc. #:	SOP340543	Revision:	A	Effective Date:
				03/29/2019

5.0 ROLES AND RESPONSIBILITIES

Laboratory Director/Supervisor The Laboratory Director/Supervisor directs laboratory operations, supervises technical personnel and reporting of findings, and is responsible for the proper performance of all laboratory procedures. The Laboratory Director/Supervisor oversees the personnel who follow the SOPs within the laboratory and is responsible for ensuring the personnel are certified and have sufficient experience to handle clinical samples.

Certified Assay Operator A Certified Assay Operator may be a Laboratory Technician/Technologist, Research Associate, or Laboratory Scientist who has been certified through DCTD training on this SOP. The Certified Assay Operator works under the guidance of the Laboratory Director/Supervisor. This person performs laboratory procedures and examinations in accordance with the current SOP(s), as well as any other procedures conducted by a laboratory, including maintaining equipment and records and performing quality assurance activities related to performance.

- 5.1 It is the responsibility of the Laboratory Director/Supervisor to ensure that all personnel have documented training and qualification on this SOP prior to the actual handling and processing of samples from clinical trial patients. The Laboratory Director/Supervisor is responsible for ensuring the Certified Assay Operator running the SOP has sufficient experience to handle and analyze clinical samples.
- 5.2 The Certified Assay Operator for this SOP should be well versed and comfortable with the operation of the Bond-RX™ System.
- 5.3 Digital versions of the Bond-RX™ slide information and staining process should be printed, including the slide event log and the first page of the slide detail log. The printed logs must be attached to the Batch Record in order to maintain a complete audit trail.
- 5.4 The Certified Assay Operator responsible for conducting the assay is to follow this SOP with associated addendum and complete the required tasks and associated documentation. The Batch Record ([Appendix 1](#)) must be completed in *real-time* for each experimental run, with each page *dated and initialed*.
- 5.5 All responsible personnel are to check the DCTD Biomarkers website (<https://dctd.cancer.gov/ResearchResources/ResearchResources-biomarkers.htm>) to verify that the most recent version of this SOP is being used.

DCTD Standard Operating Procedures (SOP)

Title:	γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides			Page 6 of 29
Doc. #:	SOP340543	Revision:	A	Effective Date:
				03/29/2019

6.0 MATERIALS AND EQUIPMENT REQUIRED

- 6.1** PADIS/IQC-Supplied Critical Reagents and Calibrator/control slides
- 6.1.1 Anti-phospho-histone H2AX (Ser139), clone JBW301, biotin conjugate (γH2AX-Biotin); (PADIS/IQC, Part#: 30000)
 - 6.1.2 Anti-phospho-NBS1 (Ser343), clone EP178, Digoxigenin (DIG) conjugate (pNBS1-DIG); (PADIS/IQC, Part#: 30020)
 - 6.1.3 Anti-Betacatenin Alexa Fluor® 546: Anti-Betacatenin (CTNNB1), rabbit Ab clone E247 (Abcam) custom conjugated to Alexa Fluor 546 (AF546) (Life Technologies); (β-Cat-AF546); (PADIS/IQC, Part#: 30025)
 - 6.1.4 Streptavidin, Alexa Fluor 488 conjugate; (SA-AF488); (PADIS/IQC, Part#: 40000)
 - 6.1.5 Alexa Fluor 647 IgG Fraction Monoclonal Mouse Anti-Digoxin; (DIG-AF647); (PADIS/IQC, Part#: 40008)
 - 6.1.6 γH2AX Cal/Con Slides; (PADIS/IQC, Part#: 60000)
 - 6.1.7 pNBS1 Cal/Con slides; (PADIS/IQC, Part#: 60006)
- 6.2** DAPI dihydrochloride, FluoroPure™ grade (Invitrogen, Cat#: D21490)
- 6.3** Pipettes (100-1000 μL, 50-200 μL, 2-20 μL, 0.2-2 μL) and tips
- 6.4** 50-mL polypropylene tubes (e.g., Becton Dickinson, Cat#: 352098)
- 6.5** Premium cover glasses, approx. 50 mm x 22 mm (e.g., Fisher Scientific, Cat#: 12-548-5E; Thermo Scientific; Cat#: 12440S)
- 6.6** Kimwipes (e.g., Fischer Scientific, Cat#: 06-666A)
- 6.7** Slide mailer/folder (e.g., Leica Microsystems, Cat#: 3802617)
- 6.8** Sterile-filtered, molecular biology grade deionized (DI) water (e.g., Invitrogen, Cat#: 10977-015)
- 6.9** 10X phosphate-buffered saline (PBS; e.g., Invitrogen, Cat#: 70013-073) [Dilute 1:10 in DI water to prepare 1X PBS for use in assay.]
- 6.10** Anhydrous ethanol, histology grade (Fisher Scientific, Cat#: A405-20 [Filtered using 0.22 μm pore size before use.]) ACS/USP Grade can be purchased and used without filtration (Pharmco-AAPER, Cat#: 111000200PL05)
- 6.11** Xylene, ACS grade (e.g. EMD Millipore, Cat# XX0055-3)
- 6.12** ProLong® Gold antifade reagent (Invitrogen, Cat#: P36930)
- 6.13** Bond-RX™ Autostainer (Leica Microsystems, Cat#: 21.2701)
- 6.14** Bond Dewax Solution (Leica Microsystems, Cat#: AR9222)
- 6.15** Bond Epitope Retrieval Solution 2 (Leica Microsystems, Cat#: AR9640)
- 6.16** Bond Open Container – 10 pack; 30 mL (Leica Microsystems, Cat#: OP309700); alternate container sizes are listed in Batch Records.
- 6.17** Bond Research Detection Kit (Leica Microsystems, Cat#: DS9455)
- 6.18** Bond Primary Antibody Diluent (Leica Microsystems, Cat#: AR9352)
- 6.19** Normal Goat Serum Blocking Solution (e.g. Vector Laboratories, Cat#: S-1000)
- 6.20** Bond Universal Covertiles, 160 Pack (Leica Microsystems, Cat#: S21.4611)
- 6.21** Bond Wash Solution 10X Concentrate (Leica Microsystems, Cat#: AR9590)
- 6.22** Bond Universal Slide Labels and Printing Ribbon kit (Leica Microsystems, Cat#: S21.4564.A)

DCTD Standard Operating Procedures (SOP)

Title:	γ H2AX, pNBS1 IFA Staining with β -Catenin Segmentation for Tumor Biopsy Slides			Page 7 of 29	
Doc. #:	SOP340543	Revision:	A	Effective Date:	03/29/2019

- 6.23 Tissue-Tek[®] Slide Staining Dish White (Sakura, Cat#: 4457)
- 6.24 Tissue-Tek 24-Slide Holder with Detachable Handle (Sakura, Cat#: 4465)
- 6.25 -80°C and -20°C freezers
- 6.26 2°C -8°C refrigerator
- 6.27 Clinical slides prepared following SOP340550 with paraffin-embedded biopsy samples and control tissues on each slide

DCTD Standard Operating Procedures (SOP)

Title:	γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides			Page 8 of 29
Doc. #:	SOP340543	Revision:	A	Effective Date: 03/29/2019

7.0 OPERATING PROCEDURES

- 7.1** Record the name of the Certified Assay Operator, the facility running the SOP and the serial number of NCI property tag number of the Bond RX unit being utilized for this run in the Batch Record ([Appendix 1](#)).
- 7.2** If slides were shipped from a separate site, save the clinical shipping manifest for the laboratory record and attach a copy to the Batch Record.
- 7.3** Record the unique Patient ID(s) and CTEP/Protocol #(s) for the slides being assayed on each page of the Batch Record. This SOP and its associated Batch Record are sufficient for up to three Bond-RX slide trays containing up to 30 slides. A minimum of one γH2AX and one pNBS1 control slide should be included with each clinical slide run. *(If more than one set of patient slides were stained during a single run on the Bond-RX™, all patient ID(s) & protocol #(s) must appear on all pages of the batch record. If necessary, individual patient batch records can be generated by making copies of the original Batch Record.)*
- 7.4** Prior to beginning the assay, read this SOP and ensure sufficient materials and reagents are in stock to perform the assay. All reagents are to be prepared for use in one experimental run, and only in the amounts required for the specific assay.

7.5 Critical Reagents

- 7.5.1** Record the date of receipt, lot numbers, stock/supplied reagent concentration, recommended working concentration, recommended dilution and expiration dates for the critical reagents in the Batch Record ([Appendix 1, Section 1A](#)).
- 7.5.2** All Critical Reagents are to be labeled with the date of receipt and stored under the specified conditions for no longer than the recommended duration.
- Storage conditions and retest dates for all Critical Reagents are provided on the shipping manifest that accompanies the critical reagent shipment.
 - If the critical reagents are purchased directly from the manufacturer, Certified Assay Sites must qualify the reagents prior to use in the Assay. Lot-to-lot differences, particularly for primary antibodies, are expected.
- 7.5.3** **Anti-γH2AX (γH2AX-Biotin):** Mouse monoclonal antibody Clone JBW301 conjugated to biotin and supplied in PBS with sodium azide and glycerol. Concentration of the material as well as the recommended working concentration will be provided by lot.
- 7.5.4** **Anti-NBS1 pS343 Digoxigenin Conjugate (pNBS1-DIG):** Rabbit monoclonal antibody Clone EP178 (Abcam) custom conjugated to DIG supplied in PBS with sodium azide and glycerol. Concentration of the material as well as the recommended working concentration will be provided by lot.
- 7.5.5** **Anti-βcatenin-Alexa Fluor 546 (β-Cat-AF546):** Rabbit monoclonal clone E247 (Abcam) custom conjugated to Alexa Fluor 546 (AF546) is supplied as a stock solution in PBS with sodium azide. The concentration of the material as well as the recommended working concentration of the material will be provided by lot.
- 7.5.6** **Streptavidin Alexa Fluor 488 (SA-AF488):** Conjugate supplied as a stock solution in PBS with BSA, sodium azide, and glycerol. The concentration of the material will be provided by lot. The optimal working concentration of SA-AF488 is 10 μg/mL.

DCTD Standard Operating Procedures (SOP)

Title:	γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides			Page 9 of 29
Doc. #:	SOP340543	Revision:	A	Effective Date:
				03/29/2019

- 7.5.7 **Alexa Fluor 647 IgG Fraction Monoclonal Mouse Anti-Digoxin (DIG-AF647):** Conjugate supplied as a stock solution in PBS with BSA, sodium azide, and glycerol. Concentration of the material will be provided by lot. The optimal working concentration of DIG-AF647 is 17 μg/mL.
- 7.5.8 **Control Slides** should be stored in a desiccator at 2°C to 8°C away from volatile chemicals. Each clinical staining run should include a minimum of one γH2AX and one pNBS1 control slide.
- 7.5.9 **DAPI stock solution** is a 14.3 mM (5 mg/mL) solution in DI water. Aliquots can be stored at -20°C for up to 1 year. The thawed aliquots can be stored at 4°C for up to 3 months. The solution is light sensitive, so it should be protected from light.
- 7.6 If not already done, program the following information into the Bond-RX™ System prior to experimental setup:
- 7.6.1 Facility or laboratory running the assay should be added to the “Researchers List”.
- 7.6.2 To enter any new antibodies and Bond Open Containers see [Appendix 2A & 2B](#) for instructions.
- 7.6.3 To enter or to verify the staining protocol, see [Appendix 2, Section 3](#).
- 7.6.4 Verify that the HIER Protocol “HIER 10 min with ER 2” matches that listed in [Appendix 2, Section 3C](#).
- 7.6.5 If a new Research Detection Kit is being used, scan the bar code to open the **Add Reagent** dialog box. Select the name of the reagent from the **Reagent name** drop-down list (select “Wash Buffer” for the Open Container) and in the expiration selection put a future date (suggested 1 year from current date).
- 7.7 **Control and Clinical Slides**
- 7.7.1 For clinical slide runs, one γH2AX and one pNBS1 control slide is recommended for each Bond-RX™ slide tray used.
- 7.7.2 Clinical samples for this assay will be frozen needle biopsies collected according to SOP340507 and formalin-fixed, paraffin-embedded and sectioned according to SOP340550. A minimum of two slides must be analyzed in order to report biomarker data, and normally three to four slides are recommended for staining and analysis for each patient slide set to ensure an adequate nuclei count is achieved for each clinical biopsy. When possible, the slides are positioned in the slide trays so that a single patient slides are contained within one Bond-RX slide tray.
- 7.8 **Preparation of Reagents**
- 7.8.1 During reagent preparation, record the lot number/serial number, expiration date, and preparation date in the Batch Record ([Appendix 1, Section 1B](#)). All reagents should be labeled with date of receipt and stored under the specified conditions for no longer than the recommended durations.

Note: Some of the following reagents may be prepared ahead of time.

DCTD Standard Operating Procedures (SOP)

Title:	γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides			Page 10 of 29
Doc. #:	SOP340543	Revision:	A	Effective Date:
				03/29/2019

7.8.2 1X Bond Wash Solution

7.8.2.1 Make 1 L of 1X solution by adding 100 mL Bond 10X Wash Solution to 900 mL of DI water. Mix the solution until it is homogenous and label the bottle as “1X Bond Wash Solution” with the lot number and preparation date. Store Bond 1X and 10X Wash Solutions at 2 °C to 8°C out of direct sunlight. 1X Bond Wash Solution can be used for 4 months.

7.8.2.2 When ready for use, 1X Bond Wash Solution can be poured into the bulk container marked “Wash Buffer” located within the Bond-RX Processing Module.

7.8.3 Research Detection Kit

7.8.3.1 Add 30 mL of 1X Bond Wash Solution to the 30-mL Open Container in the kit. Note: This container is required to be loaded with the Research Detection Kit, and only used for the first Bond wash.

7.8.4 Make sure that all required bulk reagent containers have sufficient volumes before starting the Bond-RX staining procedure. The bulk reagents containers should be at least a quarter full.

7.8.4.1 The bulk reagents include: 1X Bond Wash Solution, Bond Dewax solution (only needed if on-line dewax will be utilized), anhydrous ethanol, DI water, Bond Epitope Retrieval (ER) Solution 2.

7.8.4.2 When not in use, 1X Bond Wash Solution and ER Solution containers are stored in a 2°C to 8°C refrigerator, and the other bulk reagent containers are stored in the Bond-RX™ bulk reagent cavity.

(Pre-warming the solutions that were stored in the refrigerator is not required; temperature does not adversely affect staining.)

7.8.5 **Visually inspect all solutions to be used for the assay** to ensure there is no cloudiness or precipitate present. If they are cloudy or have a precipitate, discard the solutions and clean the bottles with a mild bleach solution. Rinse the containers thoroughly with water before reuse.

7.9 **Preparation of Antibody and Ancillary Working Solutions**

7.9.1 Label three Titration Containers or Open Containers for the assay working solutions as follows: “gH2AX_pNBS1_BCat Marker”, “DIG-AF647/SA-AF488”, and “Normal Goat Serum”. The Container size is dependent on the number of slides to be stained in a run; refer to [Appendix 1, Section 3](#) for Container volumes. The Container labels correspond to the steps programmed into the staining protocol ([Appendix 2, Section 2](#)).

7.9.2 Record the lot number and expiration date of the Bond Primary Antibody Diluent and the Bond Wash in the Batch Record ([Appendix 1, Section 1B](#)).

7.9.3 Perform the calculations in the Batch Record ([Appendix 1, Section 3](#)) to prepare the working solutions as follows:

7.9.3.1 γH2AX-Biotin /pNBS1-DIG/β-Cat-AF546 Antibody Cocktail Working Solution

DCTD Standard Operating Procedures (SOP)

Title:	γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides			Page 11 of 29
Doc. #:	SOP340543	Revision:	A	Effective Date:
				03/29/2019

- The γH2AX-Biotin /pNBS1-DIG/β-Cat-AF546_Ab Cocktail should be prepared fresh using Bond Primary Antibody Diluent. This will be used as “Marker” in the staining protocol.
- The antibody stock concentration and recommended dilution for each of the three antibodies will be provided by lot.
- To ensure there is sufficient volume for all of the slides to be stained, perform the calculations in the Batch Record ([Appendix 1, Section 3](#)).
- Briefly warm each supplied antibody Critical Reagent vial and then pipette the calculated volumes of each antibody and Bond Primary Antibody Diluent into the “Marker” Container; record the preparation date in the Batch Record ([Appendix 1, Section 3](#)).

7.9.3.2 DIG-AF647/SA-AF488 Antibody Working Solution

- **NOTE: Although typical staining procedures use a single application of secondary reagents, this secondary reagent working solution has been shown to perform better with two applications.**
- The DIG-AF647/SA-AF488 Working Solution should be prepared fresh in 1X Bond Wash. This will be used as "DIG-AF647/SA-AF488" in the staining protocol.
- The antibody stock concentration for the DIG-AF647 Antibody and the SA-AF488 conjugate will be provided by lot. The optimal working concentration of SA-AF488 is 10 μg/mL and the optimal working concentration of DIG-AF647 is 17 μg/mL.
- To be sure there is sufficient volume for all of the slides to be stained (with two applications as noted), perform the calculations in the Batch Record ([Appendix 1, Section 3](#)).
- Briefly warm the vial of DIG-AF647 antibody and SA-AF488 conjugate, and then pipette the calculated volumes of each and 1X Bond Wash into the “DIG-AF647/SA-AF488” Container; record the preparation date in the Batch Record ([Appendix 1, Section 3](#)).

7.9.3.3 Normal Goat Serum Block Working Solution

- The Normal Goat Serum Working Solution should be prepared fresh to a 2% final concentration in 1X Bond Wash Solution.
- Briefly warm the Normal Goat Serum Block vial, and then pipette the calculated volumes of Normal Goat Serum and 1X Bond Wash Solution into the “Normal Goat Serum” Open Container; record the preparation date in the Batch Record ([Appendix 1, Section 3](#)).

7.9.4 It is strongly recommended to always use fresh working solutions. Working antibody solutions can be stored at 2°C to 8°C and used for up to 5 d after preparation. If you use a stored Working Solution, note this in the deviations section ([Appendix 1, Section 5](#)).

DCTD Standard Operating Procedures (SOP)

Title:	γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides			Page 12 of 29	
Doc. #:	SOP340543	Revision:	A	Effective Date:	03/29/2019

7.10 Protocol for Slide Staining in Bond-RX Processing Module

7.10.1 System Setup for Bond-RX™ Run

7.10.1.1 **Turn on** the computer and **open** the Bond software by clicking on the Bond icon, then **turn on** the Bond-RX Processing Module.

7.10.1.2 In the Bond software, select the **Slide Setup Screen**, and then select the **Add Study** button. In the **Add Study** window (as shown in the figure below), change the fields as suggested in the table below and then click **OK**.

Field	Fill in
StudyID	CTEP#-Patient ID(s) (e.g., CTEP1234-001001)**
Study Name	Date of sample processing (e.g., 2010-10-24)
Study Comments	N/A
Researcher	Facility or laboratory running assay (from drop-down list)
Dispense Volume	150 µL
Preparation Protocol	Select “ ----- “ (if off-line Dewax) (If slides were not paraffin dipped and off-line dewaxed, use *Dewax)

** A new study will need to be added for each patient so that the slide label reflects the appropriate patient.

DCTD Standard Operating Procedures (SOP)

Title:	γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides			Page 13 of 29
Doc. #:	SOP340543	Revision:	A	Effective Date: 03/29/2019

7.10.2 Add Slides to Bond-RX Run

7.10.2.1 While still in the **Slide Setup Screen**, click the **Add Slide** button, and in the **Add Slide** window (see figure below), change the fields as shown in the table below.

Slide properties ✕

2010-10-24

Study ID:
CTEP1234-001001

Researcher:

Slide ID:
0014

Study N°:
11

Study comments:

Date created:
12/17/2018 10:19:17 AM

PRE/C1D2 #5]

Tissue type: **Dispense volume:**

Test tissue 100 µL

Negative tissue 150 µL

Positive tissue

Staining mode:

Single ▼ Routine ▼

Single

Process: IHC ISH

Marker: gH2AX_pNBS1_BCat Marker ▼

Protocols

Staining: gH2AX_pNBS1_BCat Marker Panel ▼

Preparation: *---- ▼

HIER: *HIER 10 min with ER2 ▼

Enzyme: *---- ▼

OK
Copy slide
Close

DCTD Standard Operating Procedures (SOP)

Title:	γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides	Page 14 of 29
Doc. #:	SOP340543	Revision: A
		Effective Date: 03/29/2019

Field	Fill in / Select
Slide ID	Automatically generated
Study N°	Automatically generated
Study Name	Date**
Study ID	CTEP#(s) -Patient ID(s) (**)
Comments	Sample Time-point(s) & Slide # (e.g., Pre/C1D2 #5) <u>THE SLIDE NUMBER MUST BE CHANGED FOR EACH SLIDE LABEL</u>
Tissue type	Select: Test tissue (patient slides) Positive tissue (control slide)
Dispense Volume	150 uL
Staining Mode	Single Research
Process	IHC
Marker	Name of Marker(s) (e.g., gH2AX, pNBS1, β-Cat)
Staining	Staining Protocol (e.g., gH2AX_pNBS1_BCat Marker Panel)
Preparation	Select “ *---- “ (if off-line Dewax) (If slides were not paraffin dipped and off-line dewaxed, use *Dewax)
HIER	Select Epitope Retrieval Method (e.g., HIER 10 min with ER2)
Enzyme	“*----”

** Automatically generated from Add Study Screen

- 7.10.2.2 For each new slide, a **Bond Slide ID Number** will be assigned automatically and is listed in the upper left-hand corner of the window
- 7.10.2.3 For additional slides, click the **Add Slide** button at the bottom of the window. **Before adding slide change the slide #.**
- 7.10.2.4 Once all slides are entered, click **Close**.
- 7.10.2.5 Select the **Print Labels** button at the bottom of the screen to print the labels for the slides. Select **This Case** and click **OK**. If a label does not print correctly, right-click on the label and select **Print Label**.


7.10.3 Labels should appear with the following information:

Pre/D1H2 #5
CTEP8888- 34
gH2AX/pNBS1.BCat
00E7 07E
8/11/2016

- 7.10.3.1 Labels may need to be modified to get all the critical information on the label and to get it in the correct order. To modify the label, click the Bond Admin icon, select Labels, and select the PADIS label template, as shown below. Click “Activate” to use this layout then exit the Bond Admin window.

DCTD Standard Operating Procedures (SOP)

Title:	γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides			Page 15 of 29
Doc. #:	SOP340543	Revision:	A	Effective Date: 03/29/2019



7.10.3.2 Affix the printed Bond labels to the appropriate slides, and put the slides in the designated tray; make sure the labels are aligned squarely with the inside edges of the slide so that the Processing Module can scan the information.

7.10.4 Off-line Dewax of Paraffin Dipped Slides

7.10.4.1 This procedure should be followed for all Tumor Biopsy Slides and Control slides that have been dipped in paraffin to prolong stability.

- For slides that have not been dipped in paraffin, the dewaxing procedure should be carried out on the Bond-RX Automatic Staining System using the program detailed in [Appendix 2, Section 3B](#).

DCTD Standard Operating Procedures (SOP)

Title:	γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides			Page 16 of 29
Doc. #:	SOP340543	Revision:	A	Effective Date:
				03/29/2019

7.10.4.2 Prepare the reagents for the off-line dewaxing procedure in Tissue-Tek Staining dishes, and place both the clinical and Control slides in a Tissue-Tek slide rack. Deparaffinize, rehydrate, and rinse slides as follows:

Number of Containers	Volume and Reagent	Incubation Time
4	200 mL Xylene	10 min each
4	200 mL Anhydrous ethanol	3 min each
3	200 mL 95% Ethanol	3 min each
3	200 mL DI water	2 min each
1	1X Bond Wash Buffer	Final wash

7.10.4.3 Record the time of initiation of the dewaxing procedure, check off the appropriate box to acknowledge the completion of each incubation step, and record the time the slide rack is placed in the Bond Wash Buffer in the Batch Record ([Appendix 1, Section 2](#)).

7.10.5 Remove one slide from the Bond Wash and place it in the appropriate position on a Bond Slide Tray (a minimum of one γH2AX and one pNBS1 Control slide are required for each Bond-RX™ run and will normally be placed in the last position of each slide tray in use). Hold a covertile at about a 20° angle above the slide, placing the wicking end of the covertile on the bottom of the frosted end of the slide.

7.10.5.1 Using a transfer pipette gently apply 1X Bond Wash to the tip of the covertile and continue flush while carefully lowering the covertile onto the slide.

7.10.5.2 If bubbles are introduced, remove covertile, and repeat application with a fresh covertile.

7.10.5.3 Obtain next slide from the Tissue-Tek container of Bond Wash, place onto the Bond Slide Tray and repeat above covertile application process.

7.10.5.4 Load Bond Slide Trays into the Bond-RX processing modules until the trays lock.

7.10.6 Add and Load Reagents for Bond-RX™ Run

7.10.6.1 Go back to the Bond main menu and select the Reagent icon. Using the hand-held scanner, scan the Research Detection Kit and antibody working solution Containers to enter them into the Processing Module software inventory list.

- If you are using an Open Container or Research Detection Kit that is already in the Reagent list, after scanning the Container/vial the Bond-RX™ interface will report the remaining volume (inventory) in that container. If this is **sufficient volume** for your current run, proceed to the next step. Otherwise, see below.
- If the remaining volume is not sufficient, click “**Refill**” in the pop-up window before placing the containers in the Processing Module. Note: 30-mL Open Containers can only be refilled to a maximum of 40 mL volume.

DCTD Standard Operating Procedures (SOP)

Title:	γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides			Page 17 of 29
Doc. #:	SOP340543	Revision:	A	Effective Date:
				03/29/2019

- 7.10.6.2 Place the Open Containers containing the “gH2AX_pNBS1_BCat Marker”, “DIG-647/SA-AF488”, and “Normal Goat Serum” working solutions into a reagent tray, then slide the reagent tray into a reagent tray slot at the front of the machine and lock into position. These containers can also be added to the tray containing the Research Detection Kit.
- 7.10.6.3 Place the Research Detection Kit with the Open Container containing “Wash Buffer” and working solutions into the reagent tray slot at the front of the Bond-RX™ and lock into position; the Wash Buffer Container needs to be placed into the first position of the reagent tray. The Processing Module will scan the reagent container bar codes to verify loading.
- 7.10.6.4 For each loaded tray in the processing module press the **Load/Unload** button on the front of Bond RX below each tray slot to initiate scanning of the slide labels. (*Tray numbering – e.g. tray 3 will load into the Processing Module closest to the reagent trays.*)
- Note:** Once slides are loaded into the Processing Module, the staining procedure needs to be started within 15 min or new slide labels will need to be assigned.
- 7.10.6.5 Once scanned, go to the computer screen and ensure that all of the labels were read correctly. If a slide label was not read correctly, right-click the corresponding slide and manually select the **Bond Slide ID** in the window.
- 7.10.7 Once all slides and reagent containers have been scanned, the **Play** button (triangle) will activate on the **System Status Screen** on the computer. Click the **Play** button on the screen to start processing the slides. **Note:** If the **Play** button does not light up, recheck that all trays are loaded correctly and that all containers have been scanned in. An error message will be displayed on the screen. Right-click on the error message and investigate as necessary.
- Note:** If the Bond Universal Covertiles are sticking to the slides during the staining procedure (they normally slide back and forth), it is likely that there is contamination in one of the bulk reagent solutions. Discard slides and all solutions. Clean bulk reagent bottles with a mild bleach solution and then rinse thoroughly with water before reuse.

7.11 Completion of Bond-RX™ Staining Run

- 7.11.1 Allow the Prolong Gold Antifade Reagent to equilibrate to ambient temperature (*using a heat source to warm the vial is not recommended*). If the solution appears cloudy, discard according to your institution’s safety guidelines and retrieve a fresh vial. Prolong Gold should be discarded 6 months after opening.
- 7.11.2 Just prior to slide staining completion, prepare two 250 mL Tissue-Tek staining dishes.
- 7.11.2.1 Fill the first staining dish with 200 mL DI water only. Place a Tissue-Tek 24-slide holder into this diH2O staining dish.

DCTD Standard Operating Procedures (SOP)

Title:	γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides			Page 18 of 29
Doc. #:	SOP340543	Revision:	A	Effective Date:
				03/29/2019

- 7.11.2.2 In second staining dish, prepare the DAPI Working Solution by adding 10 μL of the DAPI Stock Solution to 200 mL DI water, and mix thoroughly. Protect the solution from light by covering the entire dish with aluminum foil. Record the time of Working Solution preparation in the Batch Record ([Appendix 1, Section 4](#)).
- 7.11.3 At the completion of the Bond-RX™ staining run, push the **Load/Unload** Button to unlock the slide trays; remove the trays from the Processing Module. **Note:** Once the slides are removed from Processing Module, protect from light.
- 7.11.4 One slide at a time, remove the Bond Universal Covertile and immediately place the slide in the 24-slide holder immersed in the DI water staining dish.
- 7.11.5 Once all the slides are immersed in the DI water containing staining dish, transfer the rack to the DAPI Working Solution staining dish.
- 7.11.6 Incubate the slides for 50 min at ambient temperature in the dark (cover entire dish with aluminum foil) and gently agitate every 15 min. Record the DAPI staining start time in the Batch Record ([Appendix 1, Section 4](#)).
- 7.11.7 During the incubation time, fill three additional 250-mL staining dishes with 200 mL DI water each.
- 7.11.8 After the 50 min DAPI incubation step, remove the slide rack from the DAPI Working Solution and place it into a staining dish containing fresh DI water for 5 min. Record the time slides are removed from the DAPI Working Solution in the Batch Record ([Appendix 1, Section 4](#)).
- 7.11.8.1 Repeat the DI water wash process two additional times using a fresh DI water staining dish. Confirm the completion of wash steps in the Batch Record ([Appendix 1, Section 4](#)).
- 7.11.9 One slide at a time:
- 7.11.9.1 Transfer the slides to a paper towel, and use a Kimwipe to wick away any residual liquid, taking care not to touch the tissue or let it dry out.
- 7.11.9.2 Using a 1000 μL pipette, place no more than two drops of Prolong Gold Antifade Reagent onto the sections and cover with a cover slip.
- 7.11.10 Place the slides in a slide book, lying flat in a safe location. Allow the slides to cure overnight in the dark at ambient temperature.
- 7.12** Slides should be stored in the dark at 2°C to 8°C and imaged **within 72 h** after cover slipping.
- 7.13** Review and finalize the Batch Record and document **ANY** and **ALL** deviations from this SOP during the slide staining process in the Batch Record ([Appendix 1, Section 5](#)).
- 7.14** The Laboratory Director/Supervisor should review the Batch Record and sample reports and sign the Batch Record affirming the data contained within the reports are correct ([Appendix 1, Section 7](#)).

DCTD Standard Operating Procedures (SOP)

Title:	γ H2AX, pNBS1 IFA Staining with β -Catenin Segmentation for Tumor Biopsy Slides			Page 19 of 29	
Doc. #:	SOP340543	Revision:	A	Effective Date:	03/29/2019

7.15 Clean-up

- 7.15.1 If this is the last experimental run of the day, be sure to **turn off** the Bond-RX™ Processing Module; this will ensure the lines are cleaned at the beginning of each new day when the module is turned back on. Empty the waste containers as needed.
- 7.15.2 Store ER Solutions and 1X Wash Solution bulk reagent bottles at 2°C to 8°C. The rest of the bulk reagent containers can remain inside the body of the Bond-RX Processing Module.
- 7.15.3 Bond Open Containers can be rinsed and used 3 times (90 mL total) for the **same** reagent. It is recommended to always use fresh working solutions, but working antibody solutions can be stored at 2°C to 8°C and used for up to 5 d after preparation.
- 7.15.4 Place the Bond Universal Covertiles into anhydrous ethanol for 10 min to clean. Remove from ethanol and dry with a Kimwipe for reuse. If cracked or damaged, discard.
- 7.15.5 Make sure all Bond-RX™ daily maintenance procedures have been completed. For overall maintenance, clean the bulk reagent bottles with a mild bleach solution every 3-6 months; rinse thoroughly with water before reuse. Additionally, at least once per month perform Cleaning and Maintenance as outlined in the Leica Bond-RX™ User Manual.

DCTD Standard Operating Procedures (SOP)

Title:	γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides			Page 20 of 29	
Doc. #:	SOP340543	Revision:	A	Effective Date:	03/29/2019

8.0 OPTIONAL: SHIP TO CERTIFIED ASSAY SITE FOR ANALYSIS

If the IFA analysis will be performed at a separate certified assay site, ship the slides as as stated below.

IMPORTANT: Include a copy of the Batch Record for all samples being shipped, together with the Shipping Manifest.

- 8.1 Send an e-mail to the certified assay site prior to shipping to advise recipient of scheduled shipping time. Be sure to request and receive a confirmation e-mail prior to shipping.
- 8.2 Generate a shipping list containing all the specimen records using the Shipping Manifest template as shown in [Appendix 3](#).
- 8.3 **Verify** that the contents of the package match the Shipping Manifest.
- 8.4 Print and attach the shipping address to outside of the shipping container.
- 8.5 Record the shipping date, time, tracking number, and shipping information in the Batch Record ([Appendix 1, Section 6](#)).
- 8.6 Ship the specimens **with a copy of** the Shipping Manifest and copies of the completed Batch Records for all patient specimens. Retain copies of the completed Shipping Manifest and Batch Records in your records.

DCTD Standard Operating Procedures (SOP)

Title:	γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides			Page 21 of 29
Doc. #:	SOP340543	Revision:	A	Effective Date:
				03/29/2019

APPENDIX 1: BATCH RECORD

NOTE: Record times using **military time** (24-h designation); for example, specify 16:15 to indicate 4:15 PM

Certified Assay Operator: _____

Facility/Laboratory Running Assay: _____

Serial or NCI Property Tag Number of Bond RX: _____

Patient ID	CTEP/ Protocol #	Slide #

1. Reagents

A. Critical Reagents

Reagent Name	Date Received/ Prepared	Lot Number	Stock Reagent Conc.	Recommended Working Conc.	Recommended Dilution	Retest Date
Anti-γH2AX-Biotin, Part: 30000	/ /					/ /
Anti-pNBS1-DIG, Part: 30020	/ /					/ /
Anti-β-Catenin-AF546, Part: 30025	/ /					/ /
Anti-DIG-AF647, Part: 40008	/ /			17 μg/mL		/ /
SA-AF488, Part: 40000	/ /			10 μg/mL		/ /
DAPI	/ /		14.3 mM (5 mg/mL)	0.25 μg/mL	1: 20,000	/ /
gH2AX Cal/Con Slides, Part: 60000	/ /			N/A	N/A	/ /
pNBS1 Cal/Con Slides, Part: 60006	/ /			N/A	N/A	/ /

BATCH RECORD: INITIALS _____ DATE: _____

DCTD Standard Operating Procedures (SOP)

Title:	γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides			Page 22 of 29
Doc. #:	SOP340543	Revision:	A	Effective Date: 03/29/2019

Patient ID(s):

CTEP/Protocol ID(s):

B. Reagent Log

Reagent	Stock Solution		Working Solution	
	Lot#	Expiration Date	Concentration	Preparation Date
10X Bond Wash Solution		/ /	1X Solution	/ /
Bond Primary Antibody Diluent		/ /	N/A	N/A

2. Off-line Dewaxing of Paraffin Dipped Slides

Were slides dipped in paraffin to prolong stability? Yes No

A. Off-line Dewax Reagent Applications

Record the times and acknowledge the reagent applications step below:

Step	Time
Time Off-line Dewax Procedure Began	:
Four, 10 min Xylene Incubations Completed	<input type="checkbox"/>
Four, 3 min Anhydrous Ethanol Incubations Completed	<input type="checkbox"/>
Three, 3 min 95% Ethanol Incubations Completed	<input type="checkbox"/>
Three, 2 min DI H ₂ O Rinses Completed	<input type="checkbox"/>
Time Slides Placed in 1X Bond Wash Solution	:

BATCH RECORD:

INITIALS _____

DATE: _____

DCTD Standard Operating Procedures (SOP)

Title:	γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides			Page 23 of 29
Doc. #:	SOP340543	Revision:	A	Effective Date: 03/29/2019

Patient ID(s):
CTEP/Protocol ID(s):

3. Preparation of Working Solutions

Reagent Containers	Max. Vol. (mL)	Minimum Required Residual Vol. (μL)
Bond Open Containers, 30 mL (Leica Microsystems, Cat#: OP309700)	30	1500
Bond Titration Kit (Containers and Inserts; Leica Microsystems, Cat#: OPT9049)	6	300
Bond Open Containers, 7 mL (Leica Microsystems, Cat#: OP79193)	7	1000

Reagent Mixture Desc.	Date (Prepared)	Reagent Name	Diluent Name	A. Suggested Dilution	B. Total # <input type="text"/> slides X 300 μL	C. Residual vol (μL) 300μL, 1000μL, or 1500μL (see minimum requirement in the table above and enter below)	D. Total Volume: (needed for staining) B+C (μL)	E. Vol. reagent: total volume/ dilution factor D ÷ A (μL)	F. Vol. of Diluent D – M (or R)* (μL)
Primary Antibody Mixture		1. Anti-γH2AX Ab Marker 1	Bond Primary Antibody Diluent						*(Marker Ab) D – (M1+M2+M3)
		2. Anti-pNBS1-DIG Ab Marker 2							
		3. Anti-Betacatenin-AF546 Ab Marker 3							
Secondary Reagent Mixture		Anti-DIG-647 Reagent 1	1X Bond Wash						*Secondary Reagent D – (R1+R2)
		SA-AF488 Reagent 2							

BATCH RECORD: INITIALS _____ DATE: _____

DCTD Standard Operating Procedures (SOP)

Title:	γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides			Page 24 of 29
Doc. #:	SOP340543	Revision:	A	Effective Date: 03/29/2019

Patient ID:

CTEP/Protocol ID:

Date (Prepared)	Ancillary Reagent Name	Diluent Name	A. Suggested Dilution	B. Total # <input type="text"/> slides X 150 μL	C. Residual vol (μL) – 300μL, 1000μL, or 1500μL (see minimum requirement in the table above and enter below)	D. Total Volume: (needed for staining) B+C (μL)	E. Vol. reagent: total volume/ dilution factor D ÷ A (μL)	F. Vol. of Diluent D – E (μL)
	Goat Serum Block	1X Bond Wash	1:50					

After Bond-RX™ staining run is complete, print the Run Event log and the first page of the Run Detail Log*. Attach the documents to the Batch Record. *if there was an adverse event during the Bond-RX™ staining run, the entire Run Detail Log should be printed and attached to the Batch record.

4. Staining of Slides

A. DAPI Staining and Cover Slip Application

Just prior to staining with DAPI, prepare DAPI Working Solution by diluting 10 μL DAPI stock (5 mg/mL) with 200 mL DI water in a 250-mL staining dish. Discard excess Working Solution at end of the assay run.

	Time
Slide Trays Removed from Processing Module	:
DAPI Working Solution Added to Slides	:
DAPI Working Solution Removed	:
Three, 5 min DI Water Washes Completed	<input type="checkbox"/>
ProLong Gold Antifade Reagent with Cover Slips Added	:

BATCH RECORD:

INITIALS _____

DATE: _____

DCTD Standard Operating Procedures (SOP)

Title:	γ H2AX, pNBS1 IFA Staining with β -Catenin Segmentation for Tumor Biopsy Slides			Page 25 of 29
Doc. #:	SOP340543	Revision:	A	Effective Date: 03/29/2019

Patient ID:

CTEP/Protocol ID:

5. Notes, including any deviations from the SOP:

6. Shipping to Certified Assay Site

Date and time samples shipped: _____

Tracking information: _____

Attach copy of Shipping Manifest

7. Laboratory Director/Supervisor Review of Batch Record

Laboratory Director/Supervisor: _____ (PRINT)

_____ (SIGN)

Date: _____

BATCH RECORD: INITIALS _____ DATE: _____

DCTD Standard Operating Procedures (SOP)

Title:	γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides			Page 26 of 29
Doc. #:	SOP340543	Revision:	A	Effective Date:
				03/29/2019

APPENDIX 2: BOND-RX PROCESSING MODULE

1. Modifications to SOP for running a single slide tray in Bond-Max System

When using a Bond Titration Container with Insert, scan the bar code on the titration container when programming the Bond-RX System, and clearly label each container as Marker or DIG-AF647. The Bond Container Insert should be discarded after use, but the Bond Titration Container can be reused multiple times.

2. Register new antibodies and Open Containers in the Bond-Max System

- A. On the Reagent Screen, add “gH2AX_pNBS1_BCat Marker”, “DIG-AF647/SA-AF488” and “Goat Serum Block” to the reagent list as follows:

Field	gH2AX_pNBS1_BCat Panel	Goat Serum Block	DIG-AF647/SA-AF488
Name:	<i>gH2AX_pNBS1_BCat Marker</i>	<i>Normal Goat Serum</i>	<i>DIG-AF647/SA-AF488</i>
Abbreviated name:	<i>DDR Triplex</i>	<i>NGS</i>	<i>DIG-AF647/SA-AF488</i>
Type:	<i>Primary</i>	<i>Ancillary</i>	<i>Ancillary</i>
Single/double stain	<i>Double</i>	<i>Double</i>	<i>Double</i>
Default Staining protocol:	<i>“gH2AX_pNBS1_BCat Panel”</i>	<i>N/A</i>	<i>N/A</i>
Default HIER protocol:	<i>HIER 10 min with ER2</i>	<i>N/A</i>	<i>N/A</i>
Default enzyme protocol:	<i>*- - -</i>	<i>N/A</i>	<i>N/A</i>
Preferred	<i>Selected</i>	<i>Selected</i>	<i>Selected</i>

- B. Scan the new Open Container, Titration Kit Container, or Research Detection Kit Container bar codes to open the Add Reagent dialog box. Select the appropriate reagent name from the Reagent name drop-down list and label the Containers with the antibody names for easy identification. Repeat this procedure for each of the assay reagents. The Containers will not need to be entered again until a new Container, and therefore new bar code, is used.

3. Staining Protocols

Create the following staining protocol (A), “gH2AX_pNBS1_BCat Marker Panel”, on the Bond-RX Processing Module. Protocols B and C are pre-programmed protocols on the Bond-RX Processing Module and will be used for the gH2AX_pNBS1_BCat Marker Panel.

DCTD Standard Operating Procedures (SOP)

Title:	γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides	Page 27 of 29
Doc. #:	SOP340543	Revision: A Effective Date: 03/29/2019

A. Staining Protocol: “gH2AX_pNBS1_BCat Marker Panel” (protocol entered by user)

Solution	Temperature	Time (min)*
Wash Buffer†	Ambient	0
Bond Wash Solution	Ambient	5
Bond Wash Solution	Ambient	0
Bond Wash Solution	Ambient	0
Normal Goat Serum	Ambient	20
gH2AX_pNBS1_BCat Marker	Ambient	30
gH2AX_pNBS1_BCat Marker	Ambient	30
Bond Wash Solution	Ambient	5
Bond Wash Solution	Ambient	5
Bond Wash Solution	Ambient	5
Bond Wash Solution	Ambient	0
DIG-647_SA-AF488	Ambient	30
DIG-647_SA-AF488	Ambient	30
Bond Wash Solution	Ambient	0
Bond Wash Solution	Ambient	5
Bond Wash Solution	Ambient	5
Bond Wash Solution	Ambient	0
Bond Wash Solution	Ambient	0
Bond Wash Solution	Ambient	0
Bond Wash Solution	Ambient	0

**A time of zero indicates that the solution is applied, but that minimal time elapses before the next application.*

† The Bond-RX Processing Module requires one established solution be used from its reagent selection list. For the Research Detection Kit, 1X Bond Wash Solution is placed into a 30-mL Open Container and is used in this protocol.

DCTD Standard Operating Procedures (SOP)

Title:	γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides			Page 28 of 29
Doc. #:	SOP340543	Revision:	A	Effective Date:
				03/29/2019

B. Preparation Protocol: “*Dewax” (using Processing Module preset protocol)

Solution	Temperature (°C)	Time
Bond Dewax Solution	72	30 sec
Bond Dewax Solution	72	0
Bond Dewax Solution	Ambient	0
100% Ethanol	Ambient	0
100% Ethanol	Ambient	0
100% Ethanol	Ambient	0
Bond Wash Solution	Ambient	0
Bond Wash Solution	Ambient	0
Bond Wash Solution	Ambient	5 min

C. HIER Protocol: “*HIER 10 min with ER 2 (using Processing Module preset protocol)

Solution	Temperature (°C)	Time (min)
Bond ER 2 Solution	Ambient	0
Bond ER 2 Solution	Ambient	0
Bond ER 2 Solution	100	10
Bond ER 2 Solution	Ambient	0
Bond Wash Solution	Ambient	0
Bond Wash Solution	Ambient	0
Bond Wash Solution	Ambient	0
Bond Wash Solution	Ambient	3

DCTD Standard Operating Procedures (SOP)

Title:	γH2AX, pNBS1 IFA Staining with β-Catenin Segmentation for Tumor Biopsy Slides			Page 29 of 29
Doc. #:	SOP340543	Revision:	A	Effective Date: 03/29/2019

APPENDIX 3: STAINED SLIDES SHIPPING MANIFEST

Ship From:		Shipping Manifest			Ship To:	
Contact Name:					Attn:	
Tel:		Tel:				
E-mail:		E-mail:				
Shipping Date:		Carrier:				
In Package	Item No.	Patient ID	Clinical Protocol/CTEP#	Bond Slide ID	Comment	
<input checked="" type="checkbox"/>	<i>Example</i>	<i>Calibrator Lot # or Patient 54</i>	<i>CTEP 8888</i>	<i>78AD</i>	<i>Calibrator slide or Patient slide (Pre & D1H2)</i>	
<input type="checkbox"/>	1					
<input type="checkbox"/>	2					
<input type="checkbox"/>	3					
<input type="checkbox"/>	4					
<input type="checkbox"/>	5					
<input type="checkbox"/>	6					
<input type="checkbox"/>	7					
<input type="checkbox"/>	8					
<input type="checkbox"/>	9					
<input type="checkbox"/>	10					