

MCCRD-SOP0064: Detection of Human Leukocyte Antigen (HLA) Typing from Whole Exome Sequencing Data

Laboratory: Molecular Characterization and Clinical Assay Development Laboratory

Revision Date: 2/3/2023

Page 1 of 4

MCCRD-SOP0064: Detection of Human Leukocyte Antigen (HLA) Typing from Whole Exome Sequencing Data

Effective Date: 2/3/2023

Please check for revision status of the SOP at

<https://pdmr.cancer.gov/sops/>

PDMR **NCI Patient-Derived Models Repository**
An NCI Precision Oncology InitiativeSM Resource

TABLE OF CONTENTS

1.0	PURPOSE/SCOPE	3
2.0	DESCRIPTION OF HLA TYPING DETECTION	3
3.0	CODE DESCRIPTION.....	3
4.0	EXAMPLE HLA PROFILE REPORTING IN THE NCI PDMR DATABASE:	4

MCCRD-SOP0064: Detection of Human Leukocyte Antigen (HLA) Typing from Whole Exome Sequencing Data

Laboratory: Molecular Characterization and Clinical Assay Development Laboratory

Revision Date: 2/3/2023

Page 2 of 4

APPROVALS

Laboratory Director:	Biswajit Das, Ph.D. Chris Karlovich, Ph.D. Mickey Williams, Ph.D.
Author/Owner:	Lily Chen, Ph.D.

VERSION INFORMATION

1. Change History

Revision	Description
	Internal SOP used by MOCHA Laboratory
2/3/2023	Standardize SOP for posting to PDMM internal site for use by designated NCI intramural laboratories

2. Related SOPs

MCCRD_SOP0011: Whole Exome Sequencing Data Analysis Pipeline and Specifications

3. Reference Documentation

Number	Title/Link
[1]	Szolek, A, Schubert, B, Mohr, C, Sturm, M, Feldhahn, M, and Kohlbacher, O (2014). OptiType: precision HLA typing from next-generation sequencing data Bioinformatics, 30(23):3310-6.
[2]	https://github.com/FRED-2/OptiType

MCCRD-SOP0064: Detection of Human Leukocyte Antigen (HLA) Typing from Whole Exome Sequencing Data

Laboratory: Molecular Characterization and Clinical Assay Development Laboratory

Revision Date: 2/3/2023

Page 3 of 4

1.0 PURPOSE/SCOPE

This Standing Operating Procedure (SOP) describes procedures for detection of Human Leukocyte Antigen (HLA) typing using whole exome sequencing (WES) data for reporting in the NCI Patient-Derived Models database as performed by the Molecular Characterization Laboratory (MoCha) at the Frederick National Laboratory for Cancer Research. **This SOP is for research purposes only and no clinical samples will be processed using this SOP.**

2.0 DESCRIPTION OF HLA TYPING DETECTION

- 2.1 The mouse reads are first filtered out from whole exome sequence (WES) data following the WES data analysis pipeline in the SOP MCCRD_SOP0011.
- 2.2 HLA reads are extracted from paired-end human only sequencing data.
- 2.3 HLA type is estimated using OptiType^[1,2].

3.0 CODE DESCRIPTION

- 3.1 HLA reads are extracted from paired-end sequencing data


```
razers3 -i 95 -m 1 -dr 0 -tc 20 -o ${sample}_1.bam
$OPTITYPE_HOME/data/hla_reference_dna.fasta ${sample}_hg19_R1.fastq.gz
razers3 -i 95 -m 1 -dr 0 -tc 20 -o ${sample}_2.bam
$OPTITYPE_HOME/data/hla_reference_dna.fasta ${sample}_hg19_R2.fastq.gz
samtools bam2fq ${sample}_1.bam >${sample}_1.fastq
samtools bam2fq ${sample}_2.bam >${sample}_2.fastq
```
- 3.2 OptiType is used to detect HLA type for PDX sample


```
OptiTypePipeline.py -i ${sample}_1.fastq ${sample}_2.fastq --dna -v -o
${sample} -c $OPTITYPE_HOME/config.ini
```
- 3.3 HLA type is reported at model level with #concordance/#total PDX specimens in the model. The data for PDC and PDORG is reported separately.

MCCRD-SOP0064: Detection of Human Leukocyte Antigen (HLA) Typing from Whole Exome Sequencing Data

Laboratory: Molecular Characterization and Clinical Assay Development Laboratory

Revision Date: 2/3/2023

Page 4 of 4

4.0 EXAMPLE HLA PROFILE REPORTING IN THE NCI PDMR DATABASE:

4.1 All PDX tumors sequenced (n=5) are in concordance

HLA Profile(s)

A1	A2	B1	B2	C1	C2	PDX Concordance / PDM Type
A*02:01	A*25:01	B*08:01	B*14:02	C*07:01	C*08:02	5/5

4.2 In a small fraction of models, intra-model HLA variation in sequenced PDX tumors is observed.

HLA Profile(s)

A1	A2	B1	B2	C1	C2	PDX Concordance / PDM Type
A*01:01	A*32:01	B*08:01	B*14:01	C*05:01	C*07:01	1/6
A*01:01	A*32:01	B*08:01	B*14:01	C*07:01	C*08:02	3/6
A*01:01	A*32:01	B*08:01	B*14:01	C*07:01	C*08:04	2/6