

SOP30102: Preparation of Matrigel-Coated Flasks for Adherent Patient-Derived In Vitro Cultures

Laboratory: Patient-Derived Models Repository

Revision Date: 2/28/2022

Page 1 of 8

SOP30102: Preparation of Matrigel-Coated Flasks for Adherent Patient-Derived In Vitro Cultures

Effective Date: 2/28/2022

Please check for revision status of the SOP at

<https://pdmr.cancer.gov/sops/>

PDMR NCI Patient-Derived Models Repository
An NCI Precision Oncology InitiativeSM Resource

TABLE OF CONTENTS

1.0	PURPOSE/SCOPE	3
2.0	SAFETY	3
3.0	CLEAN-UP	3
4.0	EQUIPMENT	4
5.0	PROTOCOL 1: MATRIGEL [®] -COATED FLASKS	5
6.0	PROTOCOL 2: BME-COATED FLASKS	7

SOP30102: Preparation of Matrigel-Coated Flasks for Adherent Patient-Derived In Vitro Cultures

Laboratory: Patient-Derived Models Repository

Revision Date: 2/28/2022

Page 2 of 8

VERSION INFORMATION

1. Change History

Revision	Description
	Internal SOP used by PDMR In Vitro Laboratory
10/15/2017	Standardize SOP for posting to PDMR internal site for use by designated NCI intramural laboratories
5/14/2018	Updated reference SOPs and Purpose/Scope section
9/6/2018	Clarify steps in Matrigel coating and length of time for storage before use.
1/16/2019	Added the need for Pen/Strep in the coating solution. Streamlined protocol for readability.
2/28/2022	Added protocol for coating flasks with Basement Membrane Extract (BME)

2. Related SOPs

SOP30103: Initial Culture, Sub-culture, and Cryopreservation of Adherent Patient-Derived Tumor Cultures (PDCs)
SOP30105: Initial Culture and Sub-culture of Patient-Derived Cancer-Associated Fibroblasts (CAFs)

1.0 PURPOSE/SCOPE

This Standing Operating Procedure (SOP) describes preparation of Matrigel-coated or BME-coated plates for successful thawing and early culture of adherent Patient-Derived Tumor Cultures (PDCs) and Cancer-Associated Fibroblasts (CAFs) under BSL-2 safety criteria.

This SOP is used/performed by the Biological Testing Branch (BTB) at NCI-Frederick, Frederick National Laboratory for Cancer Research.

2.0 SAFETY

BTB treats all patient-derived in vitro cell cultures under Biosafety Level 2 (BSL2) conditions even when PCR-based screening has not detected the presence of a known set of human pathogens. All work is conducted in a biological safety cabinet (BSC) using personal protective equipment and avoiding the use of sharps where possible. All materials potentially exposed to the cell cultures are disinfected by exposure to a 10% bleach solution for a minimum of 10 minutes, double bagging for autoclaving or incineration. Consult with your facility safety professionals regarding the safe handling of BSL2 studies.

3.0 CLEAN-UP

- 3.1** All materials in contact with patient tissue, as well as the mice carrying patient tumor samples and cultures derived from patient tumor samples, are treated as a potential health threat (BSL-2 precautions) since the human tissues could retain human pathogenic agents even if they do not replicate in mouse cells (e.g., EBV, HPV, etc).
- 3.2** Flush/soak any items (e.g., tubes, syringes, petri dishes, lab mats, etc) that were in contact with human tissue with disinfectant (e.g., 10% bleach, commercial hydrogen peroxide disinfectant, 2% Virkon®) for a minimum of 10 minutes before disposal in biohazard waste or sharps containers (follow institutional guidelines and manufacturer's recommendations).
- 3.3** For items that can't be rinsed (e.g., micropipettors), wipe down thoroughly with bleach-soaked gauze or other appropriate disinfectants.

SOP30102: Preparation of Matrigel-Coated Flasks for Adherent Patient-Derived In Vitro Cultures

Laboratory: Patient-Derived Models Repository

Revision Date: 2/28/2022

Page 4 of 8

4.0 EQUIPMENT

4.1 Equipment

- 4.1.1** 50-mL, 25-mL, 10-mL, 5-mL pipettes, sterile
- 4.1.2** 15 and 50-mL polypropylene tubes, sterile
- 4.1.3** Tissue Culture flasks, sterile, vented
- 4.1.4** Pipetman and sterile tips
- 4.1.5** Waste container Bleach (Clorox, 5.25% Hypochlorite) diluted 1:10, 2% Virkon®, or similar disinfectant
- 4.1.6** Refrigerator (4°C) and freezer (−20°C)
- 4.1.7** 37°C Incubator (5% CO₂, humidified)
- 4.1.8** Biological Safety Cabinet (BSC) meeting biosafety level 2 (BSL2) standards
- 4.1.9** Personal Protective Equipment (PPE) at a minimum laboratory coat, with fitted sleeves, latex or nitrile gloves and safety glasses

5.0 PROTOCOL 1: MATRIGEL®-COATED FLASKS

5.1 Reagents

5.1.1 1X HAM's F-12 Nutrient Mix, with L-glutamine (Invitrogen, Cat#: 11765-047) with 100 U/mL Pen/Strep final (Invitrogen, Cat#: 1514022, 10000 U/mL)

5.1.2 Matrigel® Matrix

5.1.2.1 Matrigel®, High concentration (BD Biosciences, Cat#: 354248)

OR

5.1.2.2 Matrigel®, Standard concentration (BD Biosciences, Cat# 354234)

IMPORTANT: All Matrigel® purchases should be submitted specifying PCR-tested LDEV-Negative Matrigel®. If not, there is a possibility of LDEV contamination which can result in LDEV+ tumors

5.2 Prepare Matrigel® Working Solution

5.2.1 Chill pipettes and conical tubes in a -70°C freezer overnight and then place on wet ice prior to use.

5.2.2 Thaw Matrigel® overnight by placing a vial of Matrigel, buried in ice, in the refrigerator.

5.2.3 Using cold pipettes and tubes, make a Matrigel® Working Solution with 1X HAM's F-12 Nutrient Mix, with L-glutamine, + Pen/Strep.

5.2.3.1 Make a 25% Matrigel Working Solution if using High Concentration

5.2.3.2 Make a 50% Matrigel Working Solution if using Standard Concentration

5.2.4 Aliquot Matrigel® Working Solution either as 5-mL aliquots into chilled sterile 15-mL conical tubes or 1-mL aliquots into 2-mL sterile screw-capped tubes and place into a -20°C freezer.

5.3 Preparation of Matrigel®-coated flasks

5.3.1 The day before coating, remove the appropriate number of Matrigel® Working Solution aliquots (prepared in Section 5.2) from the freezer and thaw overnight buried in ice, in the refrigerator.

5.3.2 The recommended volume of 2.5% Matrigel needed per well/flask are:

Plate/Flask size	Volume 2.5% Matrigel®/well or flask
96-well plate	75 µL/well
24-well plate	0.3 mL/well
6-well plate	2 mL/well
T25 flask	3.0 mL
T75 flask	6.0 mL
T162 flask	8.0 mL

5.3.3 Prepare a 2.5% Matrigel solution using the 25% or 50% Matrigel Working Solution in 1X HAM's F-12 Nutrient Mix, with L-glutamine, + Pen/Strep. Keep solution chilled.

5.3.4 Coat Tissue Culture Plates/Flasks with 2.5% Matrigel® Solution

5.3.4.1 Matrigel® coated plates/flasks should be prepared at least 1 hour before use.

5.3.4.2 Using sterile procedures, coat the growth surface of the plate/flask using the recommended volume based on the size of the well/flask (Section 5.3.2). Be sure to rock the plate back and forth to completely coat the surface of the plate/flask.

5.3.4.3 Incubate plates/flasks for a minimum of 30 minutes at ambient temperature for polymerization. This process can be enhanced by incubating in a 37°C incubator.

5.3.4.4 Immediately before use, remove excess media from the flask/plate and discard taking care to not dislodge the Matrigel® coating.

- Flasks can be prepared up to 7 days before use and stored at 4°C following the polymerization step (leave excess media on flask/plate until use). Bring to ambient temperature or 37°C before adding cells.

6.0 PROTOCOL 2: BME-COATED FLASKS

6.1 Reagents

6.1.1 1X HAM's F-12 Nutrient Mix, with L-glutamine (Invitrogen, Cat#: 11765-047) with 100 U/mL Pen/Strep final (Invitrogen, Cat#: 1514022, 10000 U/mL)

6.1.2 Cultrex™ Basement Membrane Extract (BME), PathClear®, with Phenol Red, 8-12 mg/mL stock (R&D Systems, Cat# 3432-005-01P)

6.1.2.1 A lot-specific Certificate of Analysis is included in each Cultrex™ BME shipment noting exact protein concentration.

6.1.2.2 Phenol Red is used as a visual indicator of coated plates

6.2 Prepare Cultrex™ BME Working Solution

6.2.1 Thaw Cultrex™ BME overnight by placing a vial of BME, buried in ice, in the refrigerator.

6.2.2 Using pipettes and chilled tubes, dilute the 100X Phenol Red solution to 1X in BME. Note: pipettes can be room temperature

6.2.3 Aliquot the Cultrex™ BME solution with 1X Phenol Red either as 5-mL aliquots into chilled sterile 15-mL conical tubes or 1-mL aliquots into 2-mL sterile screw-capped tubes and place into a -70°C freezer.

6.3 Preparation of Cultrex™ BME-coated Flasks

6.3.1 The day before coating, remove the appropriate number of Cultrex™ BME solution with 1X Phenol Red aliquots (prepared in Section 6.2) from the freezer and thaw overnight buried in ice, in the refrigerator.

6.3.2 The recommended volume of 2.5% BME needed per well/flask are:

Plate/Flask size	Volume 2.5% BME/well or flask
96-well plate	75 µL/well
24-well plate	0.3 mL/well
6-well plate	2 mL/well
T25 flask	3.0 mL
T75 flask	6.0 mL
T162 flask	8.0 mL

6.3.3 Prepare a 0.46-0.48 mg/mL BME Working Solution containing ice-cold 1X Phenol Red with 1X HAM's F-12 Nutrient Mix, with L-glutamine, + Pen/Strep. Keep solution chilled.

6.3.4 Coat Tissue Culture Plates/Flasks with 0.46-0.48 mg/mL BME Working Solution

- 6.3.4.1 BME coated plates/flasks should be prepared at least 1 hour before use.
- 6.3.4.2 Using sterile procedures, coat the growth surface of the plate/flask using the recommended volume based on the size of the well/flask (Section 6.3.2). Be sure to rock the plate back and forth to completely coat the surface of the plate/flask.
- 6.3.4.3 Incubate plates/flasks for a minimum of 1 hour at ambient temperature for polymerization. This process can be enhanced by incubating in a 37°C incubator.
- 6.3.4.4 Immediately before use, remove excess media from the flask/plate and discard.
 - Flasks can be prepared up to 7 days before use and stored at 4°C following the polymerization step (leave excess media on flask/plate until use). Bring to ambient or 37°C before adding cells.