

PAR Immunoassay Frequently Asked Questions

SOP340505 – PAR IMMUNOASSAY

1. *Can assay plates be prepared in advance?*

Yes. After coating the plates in SOP Step 7.8, block the wells according to SOP Step 7.10. Plates can then be stored at 2°C to 8°C for up to 1 week. Be sure plates do not dry out.

SOP340520 – BIOPSY SPECIMEN PROCESSING

There are currently no questions for this SOP.

SOP340521 – PBMC SPECIMEN PROCESSING

1. *Why aren't protein levels used to normalize the PAR levels in PBMCs?*

Contamination of serum proteins can lead to artificially high protein values because albumin cannot be consistently removed. Therefore, for PBMCs, the cell count is used for normalization.

GENERAL

There are currently no general questions for the PAR Immunoassay.